DeClogger®

INSTRUCTIONS FOR USE

	
Importer, Importeur, Importör, Importador, Importateur, Dozovce, Maahantuja, Importatore, Έισαγωγέας, Importuotojas:	Primes Magreijplantsoen 33 - Suite 123 2595 AM The Hague The Netherlands

This manual applies to the **[REF]** numbers indicated below:

#911 #912 #913 #921 #922 #935

The following Bionix® DeCloggers are compatible with ENFit® connection: #911-Blue, #912-Yellow, #921-Brown, #922- Orange. ENFit® and NRFit™ are trademarks of GEDSA

	©2022 Bionix LLC 1670 Indian Wood Circle Maumee, Ohio 43537 United States Phone: 800.551.7096 Fax: 800.455.5678 Email: bionix@bionix.com www.bionix.com RM95-0007 2022/05/11 Rev. F
	

ENGLISH

DECLOGGER® - USER INSTRUCTIONS

Intended use: The intended use of the DeClogger is to clear clogs from enteral feeding tubes and thereby reduce the occurrence of feeding tube replacement procedures.

Indications for use: The DeClogger is indicated for use only in clearing occlusions and/or clogs in feeding tubes (G, J, and/or PEG style) in sizes from 14 Fr – 22 Fr. The DeClogger styles are indicated for use as follows:

- DG-00911: For G or J type enteral feeding tubes that are 14 – 16 French, and of a length greater than 39,5 cm (~15,5 in)
- DG-00912: For G or J type enteral feeding tubes that are 16 – 22 French, and of a length greater than 39,5 cm (~15,5 in)
- DG-00913: For G or J type enteral feeding tubes that are 20 – 22 French, and of a length greater than 39,5 cm (~15,5 in)
- DG-00921: For PEG or Replacement type enteral feeding tubes that are 14 – 16 French, and of a length greater than 21,5 cm (~8,5 in)
- DG-00922: For PEG or Replacement type enteral feeding tubes that are 18 – 24 French, and of a length greater than 21,5 cm (~8,5 in)

Protocol for use of the DeClogger®

Purpose:
To maintain or achieve patency of gastric and/or jejunostomy enteral tubes.

Policy
The DeClogger should be used every week to maintain unimpeded flow of an enteral formula. The DeClogger will also be used to achieve patency of a tube that has become clogged with semi-solid formula. The DeClogger should be disposed of after a single use.

Protocol
The following protocol should be followed according to the desired need.

INSTRUCTIONS

Blocked Jejunostomy Tube

- Determine the size of the gastric or jejunostomy tube.
- CAUTION:** Verify that the tube has not been modified, i.e. shortened.
- Select appropriate size DeClogger that corresponds to the size of the tube.
- Turn enteral feeding pump to PAUSE position and disconnect from enteral tube. (If gravity pump, clamp the tubing.)
- Insert the DeClogger to reach the blockage and slowly rotate two times in a clockwise direction then reverse and rotate two times in a counter-clockwise direction while removing it. Do not attempt to force the DeClogger through the entire blockage.
- Repeat step 4 until the stop disk of the DeClogger is reached without difficulty.
- Flush tube with 30-60 cc's of water.
- Once clear, reconnect the delivery tube and restart the enteral feeding.
- Discard DeClogger and document results.

Prophylactic Use
For prophylactic use, implement the DeClogger protocol every week depending on the patient's feeding and medication scheduling.

CAUTION
If the DeClogger becomes imbedded in the blockage and cannot be easily dislodged, the enteral feeding tube will need to be changed. DO NOT attempt to physically remove the DeClogger.

- WARNINGS:**
- If the DeClogger becomes imbedded in the blockage and cannot be easily dislodged, the enteral feeding tube will need to be changed. DO NOT attempt to physically remove the DeClogger.
 - Cross-contamination risk. Do not reuse DeCloggers as this may spread contamination from one patient to another patient.
 - The Bionix DeClogger #913-Green is NOT compatible with ENFit® feeding tube connectors.

Medical Device Reporting: Notice to Users and/or Patients in EU: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Clinical Benefits Indicated in Intended Use.

Contraindications: There are no absolute contraindications related to the use of the DeCloggers in clearing enteral feeding tube obstructions.

Residual Risk: Risk associated with the use of this product has been reduced as far as possible, but the product cannot completely eliminate potential patient or user harm arising from the following:

- Harm from mechanical hazards
- Harm from misuse, or use error
- Harm from unanticipated origins

	Medical Device		Do Not Use if Package is Damaged
	Manufacturer		Do Not Reuse Single-Use Device
	Date of manufacture		Single patient- multiple use
	Use-By Date		Consult Instructions For Use (IFU)
	Lot Code		Caution
	Reorder Number		Warning
	Serial number		Authorized representative in the European Community
	Sterilized using irradiation		European Conformity
	Do Not Sterilize		Prescription Only or "For Use by or on the order of a licensed medical professional"
	Non-sterile		Temperature limitation

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