

Lifestyle
MOBILITY AIDS

KN3000

AVAILABLE COLORS:



350 LBS
Weight Cap

OPERATING INSTRUCTIONS

This mobility device provides strength and comfort to anyone who may have an injury below the knee. It allows an even distribution of body weight by supporting half of the weight on the device and the other half on the non-injured leg. The device has a cushion to support both the injured leg and the ankle to make certain that no weight will be placed on the lower leg.

The device is designed for easy steering and stopping. Brakes are included on this.

- 1 With brake applied, placement of injured leg on knee pad should be centered side to side and positioned forward to cover the full length of the knee pad.
- 2 With injured leg on knee pad, stand as straight as possible. (Adjust the height of Knee Walker Pad as necessary.)
- 3 With injured foot pointing down, the propulsion leg should be maintained as close as possible to the knee pad.
- 4 Start out with small steps with propulsion leg while concentrating on keeping that leg next to the knee pad. With practice, you should be more comfortable with moving forward without deviating to the opposite side.



ASSEMBLY INSTRUCTIONS

- 1 Remove contents from carton.
- 2 Insert tiller into front of frame.
- 3 Raise the tiller and lock in place by sliding the lever up and then pushing to lock in place (from the knee pad looking forward).
- 4 To adjust the height of tiller slide up or down until it is in a comfortable position then lock in place with the knob and bolt.
- 5 Insert knee pad post into receptacle on frame. Set to desired height and secure in place by inserting knob and bolt on left side of frame facing forward and turning clockwise until tight.



HAND BRAKE USE

- 1 To use the brake, simply pull the lever with your fingertips toward the handle bar.
- 2 To lock the brake, pull the brake lever toward the handle bar and push down the spring loaded push button (located on top of the brake). When done correctly, push button will remain down and brake lever will be locked in place.
- 3 To release brake, simply pull lever toward the handle bar. Push button will pop up automatically.



WARNINGS

- Before using the knee walker, you should be trained by a healthcare professional or dealer.
- Engage the lock and practice bending, reaching and transfers on and off the knee walker.
- The use of the knee walker is individual to the person. You should develop your own methods for use based on your personal level of function and ability.
- Never attempt a maneuver that has not been practiced.
- Be aware of your surroundings when operating the knee walker. Look for hazards and avoid them.
- Avoid electrical cords, toys, and debris that could be a trip hazard.
- Do not use the knee walker on or near stairs or escalators.
- Use caution when moving from a carpeted area to hard surfaces or when changing levels, such as from a sidewalk to a street.
- Do not use while engaging in the use of mind-altering drugs or while suffering from dizziness.
- Use caution on uneven terrain, do not use in sand.
- Weight cap is 350 LBS.





Your Lifestyle Mobility Aids product is warranted to be free of defects in workmanship and materials for the lifetime of the original consumer.

This device was built to exacting standards and was carefully inspected prior to shipment. This warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service.

In the event of a defect covered by this warranty, we will at our option, repair or replace the device.

This warranty does not cover device failure due to misuse or negligence, or normal wear and tear. This warranty does not extend to non-durable components, such as rubber accessories, castors or grips, which are subject to normal wear and tear and need periodic replacement.

If you have any questions about your Lifestyle Mobility Aids device, please contact your authorized Lifestyle Mobility Aids dealer.



www.lifestylemobilityaids.com 