

# Turner

## Manual

etac®

10-07-09 74274H

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- EN** 1. Begin by positioning the new seat (chair, wheelchair, toilet or bed) in line with the existing seat.  
2. Slide Turner in under the user's feet.  
3. Adjust the height of the lower leg support so that its upper edge lines up just under the knee of the user.  
4. Adjust the height of the handle to a position that is comfortable for both user and carer.  
5. Grip the handle with both hands and places one foot on the pedal.  
6. Grip the handle with both hands.  
7. Encourage the user to stand up while, at the same time, counterbalancing the user's weight using a foot and hands (use your bodyweight by leaning backwards).  
8. Twist Turner so that the user is transferred to the front of the new seat.  
9. Step on the pedal, brace your weight and ask the user to sit down. Maintain eye contact with the user during the entire transfer process.

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**EN Safety and preventative maintenance:**

- Clean the rotating plate using a standard, solvent-free cleaning agent and wipe it dry.
- Turner must not be used to transport people or be tipped while someone is standing on the rotating plate.
- Keep the wheels and the space between the rotating plates free of dust, hair and suchlike. If required, wipe off any grease between the sliding surfaces and lubricate with silicone grease.
- Check that the screws and controls are tightened properly.

**Guarantee:**

5 year guarantee against material and manufacturing defects. For terms and conditions, see [www.etac.com](http://www.etac.com).

The expected life, during normal usage, is approximately 10 years.

If your chair needs reconditioning or repair, use only original parts from Etac or components of similar quality.

Etac will not be held responsible for damage or injury caused by use of non-original parts.

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Etac Sverige AB
Etac AB (export)
Etac A/S