



Homecraft All-in-One Bariatric Commode

Maximum User Weight: 700lbs/318kg

REF 7103248

The following manual provides descriptions relating to the installation of the Homecraft Bariatric All-in-One Commode.

INTENDED USE:

The Homecraft All-in-One Bariatric Commode offers versatility for all users. It can be used as a bedside commode, or over the toilet. The drop arm design allows for a smooth side transfer. The sturdy frame boasts a capacity limit of 700lbs/318kgs and is height adjustable from 19" to 23" (seat height).

INSTALLATION INSTRUCTIONS

NOTE: Before assembly, ensure that all parts and quantities match parts list below. Inspect all components for shipping damage. If any parts are determined to be missing or damaged, contact your dealer/retailer for replacements. ⚠️ **WARNING! DO NOT** use this equipment without first reading and understanding the user manual. If you are unable to understand these instructions, contact a healthcare professional, dealer, or technical person before installing. Otherwise, injury or damage may occur.

ASSEMBLY INSTRUCTIONS

1. Place the commode upside down on a flat surface. Insert each leg into the commode frame at each corner. Please ensure the commode is adjusted to the correct height and ensure all pins have locked in place (you will hear a clicking noise) once the pins have fully protruded through the holes.
2. Insert back into the slots, which are at the rear of the seat, and insert back until fully engaged and the pins are fully protruded.
3. The commode comes with a potty and lid for use as a commode or a splash guard for use over the toilet. Select and place on the underneath rim to secure in position.
4. To secure the armrests in position, rotate until the armrests are at the side of the commode and the armrest pin is locked into place (you will hear a clicking noise) once the pins have fully protruded through the holes.



CAUTIONS

- Ensure all 4 legs are at the correct height.
- Ensure the seat is secure before applying any weight.
- Only use the commode for its intended use. They are not replacements for ladders and are not exercise equipment or toys.
- Under no circumstances should you stand on the commode
- The commode is not intended for use during examinations or treatments.
- Use only spare parts and accessories released by Performance Health.
- Do not carry out structural changes which are not approved by the manufacturer.
- The commode must not be used in environments in which flammable gases or vapors (e.g. from anesthetics) can be expected.

CLEANING

- Small stains and dust are best removed with moist wipes.
- Larger stains can be removed by non-scouring, mild cleaning agents.
- Do not use tools with sharp edges (knives, metal scrapers) or aggressive solvents for cleaning.

MAINTENANCE

The bariatric commode is maintenance-free, however, should be checked daily basis for damage or looseness.

Carry out the following test

- Visual check of all parts for deformation, damage, wear and tear and corrosion.

The commode must not be used further if faults/defects are detected during use.

DISPOSAL

The bariatric commode consists of metal and plastic parts. Disposal of these materials must be carried out according to the regional environmental and disposal provisions.



RESTRICTIONS OF USE!

- Patients weighing more than 700lbs must not use the Homecraft bariatric all-in-one commode

LIMITED PRODUCT WARRANTY

Performance Health Supply, LLC warrants that all products shall be free from manufacturing defects for a period of one year from the date of shipment. OTHER THAN THE FOREGOING WARRANTY, PERFORMANCE HEALTH SUPPLY, LLC MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, AS TO MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER MATTER RELATING TO THE PRODUCT. Customer must notify Performance Health Supply, LLC of any manufacturing defect within five days of discovery of the defect. Upon a determination by Performance Health Supply, LLC that the product is defective, Performance Health Supply, LLC sole and exclusive obligation shall be to provide customer with a replacement product or a full refund of the purchase price, at the sole discretion of Performance Health Supply, LLC, and to reimburse customer for any return freight costs. In no event shall Performance Health Supply, LLC be liable for any other costs or expenses including lost profits or direct, indirect, consequential, incidental or punitive damages. The warranty set forth herein shall govern all purchases, notwithstanding any inconsistent terms and conditions of customer's purchase order.

Symbols Glossary located at <http://www.performancehealth.com/medical-device-glossary>

Distributed in the USA by:
Performance Health Supply, LLC
W68 N 158 Evergreen Blvd.
Cedarburg, WI 53012 USA
Phone: 800-323-5547
www.performancehealth.com

Distributed in Canada By:
Performance Health Canada Inc
6715 Millcreek Drive, Unit 2
Mississauga, ON L5N 5V2 Canada
Phone: +1800-665-9200
www.performancehealth.ca



GENTRY WAY CO., LTD.

No. 152, WuFu Road, YanShui DIST.,
Taiwan City, Taiwan

