ASP™

REF 22591

CIDEZYME™ XTRA

Multi Enzymatic Detergent for Automated Systems and Manual Cleaning

en

CIDEZYME™ XTRA

Multi Enzymatic Detergent for Automated Systems and Manual Cleaning

INSTRUCTIONS FOR USE

CIDEZYME XTRA Multi Enzymatic Detergent for Automated Systems and Manual Cleaning when used as directed is a low-foaming, mild pH cleaning solution which combines chemical and multi-enzymatic actions to rapidly remove medical soils (e.g. blood). The product also contains corrosion inhibitors to ensure compatibility with medical devices. CIDEZYME XTRA Multi Enzymatic Detergent is uniquely formulated for use in automatic cleaning systems (e.g. endoscope cleaner and reproprocessors, automatic endoscope reproprocessors) as well as for manual cleaning.

RECOMMENDED CONCENTRATIONS
For use in automatic cleaning systems, minimum 4 mL per 1 liter (1/2 oz per gallon) of water or per automatic cleaning system manufacturer's instructions.

For manual cleaning, minimum 8 mL per 1 liter (1 oz per gallon) of water.

CAUTION
As with all chemicals, good industrial hygiene practice should be followed. For further information refer to the Safety Data Sheet.

- PLEASE NOTE**
- Confirm CIDEZYME XTRA Multi Enzymatic Detergent compatibility with the instruments and equipment prior to use.
 - Do not add bleach or other chemicals.
 - Soak and clean instrument immediately after use, avoid allowing the medical soil to become dried.
 - Minimum in-use temperature for detergent is 15°C.
 - Unless the automatic cleaning system is specifically designed for elimination of manual cleaning, it will be necessary to pre-clean instruments with brushes or soft cloths as recommended by the instrument manufacturer.
 - CIDEZYME XTRA Multi Enzymatic Detergent is not intended as a solution for holding or storing equipment.

- DIRECTIONS FOR USE**
Automatic Cleaning & Disinfecting System (e.g. automatic endoscope reproprocessors or endoscope cleaner and reproprocessors)
- Connect bottle of CIDEZYME XTRA Multi Enzymatic Detergent to the machine or dispense into detergent reservoir as per system manufacturer's instructions.
 - Select a rinse cycle on an automatic system that has been validated for use with this product.
 - Ensure that the automatic rinse cycle selected will thoroughly rinse the medical device including all lumens with large volume of water equivalent to the reusable device manufacturer's recommendations.

- Manual Cleaning**
- Attach the pump to the bottle and secure tightly.
 - Add minimum 8 mL of CIDEZYME XTRA Multi Enzymatic Detergent liquid concentrate per 1 liter of water (1 oz per gallon).
 - Soak instruments for a minimum of 3 minutes. Extend soak time or use warm water for

ADVANCED STERILIZATION PRODUCTS, INC.

33 Technology Drive, Irvine, CA 92618, USA

1-888-783-7723

24-Hour Number: 1-800-451-8346 or 1-760-602-8700

equipment with dried-on organic matter. Flush solution through all channels. Soaking for long periods of time may cause discoloration of anodized aluminum.

4) Thoroughly rinse the equipment, aspirating water through all channels to remove traces of CIDEZYME XTRA Multi Enzymatic Detergent.

5) Roughly dry equipment, including channels. Equipment is now ready for reprocessing.

6) Discard used solution after each use.

STORAGE

CIDEZYME XTRA Multi Enzymatic Detergent should be stored in its original container at controlled temperature 15°- 30°C (59°- 86°F). Replace cap when not in use.

DISPOSAL

In-use concentration of CIDEZYME XTRA Multi Enzymatic Detergent can be disposed directly to drain after use or per hospital policy. After rinsing with water, the empty containers can be disposed of as non-hazardous waste or per hospital policy. Use by expiration date on the immediate container.

RISK/SAFETY INFORMATION

Danger

Contains Subtilisin, Lipase, triacylglycerol. May produce an allergic reaction. **Hazard** May cause allergy or asthma symptoms or breathing difficulties if inhaled. The mixture does not meet the criteria for classification.

Precautionary Statement

Prevention Avoid breathing mist. In case of inadequate ventilation wear respiratory protection. Observe good industrial hygiene practices.

Response

If inhaled: If breathing is difficult, remove person to fresh air and keep comfortable for breathing. If experiencing respiratory symptoms: Call a poison center/ doctor. Wash hands after handling.

Storage

Store away from incompatible materials.

Disposal

Dispose of waste and residues in accordance with local authority requirements.

**For country specific GHS and SDS information, refer to country specific SDS located at asp.com/service-and-support/sds-datasheets.

HOW SUPPLIED

Reorder

22591

Description

2 X 3.8L (1 gallon)

For a glossary of symbols used, go to eifu.asp.com.

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