

INSTRUCTIONS FOR USE

Cleaning Solutions



General Purpose Ultrasonic Cleaner Concentrate Non-Ammoniated

Instructions for use:

- Factory sealed for product reliability and protection. To open, place container on a flat steady surface and cut seal slowly. Do not squeeze container.
- **THIS PRODUCT IS A CONCENTRATE.** Dilute one part solution with ten parts water. For use directly in ultrasonic cleaning system.
- To clean, place items in accessory basket or cassette. Insert in ultrasonic tank filled with $\frac{1}{2}$ to $\frac{2}{3}$ full of properly diluted solution.
- Approximate cleaning time is 2-10 minutes depending on volume and amount of debris.
- When cleaning is complete, rinse with running water.
- Ideal as an ultrasonic carrier bath when using an auxiliary pan or beakers for other cleaning tasks.
- For optimum cleaning results, change solution at least once per day. Helps avoid cross contamination.
- **DO NOT USE ON ALUMINUM.**



Ref.#'s JEZNA

WARNING! FOR PROFESSIONAL USE ONLY. NOT FOR INTERNAL USE.

- Keep away from children.
- Use with adequate ventilation.
- Keep container closed.
- Do not mix with any other solution.
- Do not add any ingredients to this solution.
- This is a concentrate—use care to dilute properly.

Emergency and First Aid Procedures

Eye Contact: Hold eyelids open and flush with a steady, gentle stream of water for at least 15 minutes. Seek medical attention.

Skin Exposure: Wash thoroughly with water. If irritation or redness develops, seek medical attention.

Ingestion: Seek immediate medical attention. **DO NOT INDUCE VOMITING.**

Spills: Biodegradable product. Mop up spill and dispose of by dilution to a sanitary sewer system as permitted by local, state & federal regulations.

Contains: Isopropal Alcohol (CAS #68608-65-1), Potassium Hydroxide (CAS #1310-58-3)

24 Hour Emergency Telephone # INFOTRAC: 1.800.535.5053



Crosstex International, Inc.
10 Ranick Road, Hauppauge, NY 11788-4209 USA
www.crosstex.com 800.276.7783 631.582.6777

Made in USA
JEZNA-IFU Rev. 1 - 5/19

All product names are trademarks of Crosstex International, Inc., a Cantel Medical Company, its affiliates or related companies, unless otherwise noted.

CROSSTEX
A CANTEL MEDICAL COMPANY