

Product Name: **Biopsy Foam Pads**

Catalogue No.: **M476-1**

This document replaces any previous version

1. Product Description:

- **Biopsy Foam Pad:** Non-sterile, disposable rectangular foam pad cassettes

2. Packaging:

- **Case:** 10 bags of 1 000 units / 10 000 units per case

3. Product Specifications:

- **Material:**
 - Blue Polyester urethane Foam
- **Temperature range:** -40°C to +121°C
- **Autoclavable** at +121°C for up to 30 minutes

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures

5. Quality Assurance:

- No presence of contamination in foam
- Visual attributes
- Dimension measurements

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination
- Protect from any damage to the packaging

8. Recommended Use:

- Fits most standard size cassettes
- For use to hold biopsies inside histology cassettes
- check cassette's lid closure before tissue processing

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*

Product Name: Biopsy Foam Pads

Catalogue No.: M476-4

This document replaces any previous version

1. Product Description:

- **Biopsy Foam Pad:** Non-sterile, disposable rectangular foam pad for the Simport Slimsette cassettes

2. Packaging:

- **Case:** 10 bags of 1 000 units / 10 000 units per case

3. Product Specifications:

- **Material:**
 - Blue Polyurethane Foam
- **Dimensions:** 27.4 x 25.4 x 2 mm
- **Temperature range:** -40°C to +121°C
- **Autoclavable** at +121°C for up to 30 minutes

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures

5. Quality Assurance:

- No presence of contamination in foam
- Visual attributes
- Dimension measurements

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination
- Protect from any damage to the packaging

8. Recommended Use:

- Fits Simport Slimsette cassettes
- For use to hold biopsies inside histology cassettes
- check cassette's lid closure before tissue processing

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*