

TECHNICAL DATA SHEET Disinfectant Spray

Intended Use

Used to clean or disinfect hard, non-porous surface in dental offices, hospitals and other healthcare settings.

Active Ingredients

Ethyl Alcohol 72.50 % w/w

Didecyl Dimethyl Ammonium Chloride 0.33 % w/w

Other Ingredients 27.17 % w/w

Characteristic	Specification	
Material composition	Bottle	HDPE
	Child Proof Cap	HDPE
	Seal	HDPE
	Label	Film
Description	Spray bottle containing disinfectant (Refer to Appendix 1 for Specific Attribute)	

Visual



**Additional Information**

Storage Condition	Shelf Life	Country of Origin
Protect from extreme cold or heat	03 years	USA

Appendix 1 List of Item Codes and Attributes

Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Size	Color	Treatment
40066C	Medicom Pro-Surface+ Disinfectant Spray	12 Bottles/Case	32 fl. oz. (946 ml)	Clear Liquid	Non-Sterile

Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Size	Color	Treatment
40066	Medicom Pro-Surface Disinfectant Spray	12 Bottles/Case	32 fl. oz. (946 ml)	Clear Liquid	Non-Sterile