

MadaCide-FD

GERMICIDAL SOLUTION

DISINFECTANT / CLEANER / DEODORIZER

CORONAVIRUS COVID-19 UPDATE

MadaCide-FD (EPA Reg # 1130-15) has demonstrated effectiveness against **SARS-CoV-2** (the technical name for the **novel 2019 Coronavirus**) on hard, nonporous surfaces. Therefore, MadaCide-FD can be used against **SARS-CoV-2** when used in accordance with the directions for use against Rhinovirus [Type 14] on hard, non-porous surfaces. Refer to the CDC website (<https://www.cdc.gov/coronavirus/2019-ncov/index.html>) for additional information.

BROAD SPECTRUM EFFECTIVENESS

- ✓ Bactericidal
- ✓ Fungicidal
- ✓ Tuberculocidal
- ✓ Virucidal

FEATURES / BENEFITS

- ✓ Ready to Use
- ✓ Mild Scent
- ✓ Low Alcohol

ORDERING INFORMATION

No. 7020	32 oz. Spray Bottles	(12/case)
No. 7021	1 Gallon Bottles	(4/case)
No. 7022	2.5 Gallon Bottle	(each)
No. 7023	5 Gallon Bottle	(each)
No. 7024	55 Gallon Drum	(each)



No. **7020**
32 oz. Spray Bottle (12/case)

No. **7021**
1 Gallon (4/case)

MadaCide-FD

GERMICIDAL SOLUTION

MadaCide-FD is a general purpose cleaner, deodorizer and disinfectant to be used where control of cross-contamination from treated surfaces is of prime importance.

FOR USE ON ALL HARD, NON-POROUS SURFACES INCLUDING:

- Chrome
- Painted Surfaces
- Polypropylene
- Glazed Porcelain
- Formica® Brand
- Polycarbonate
- Polystyrene
- Stainless Steel
- Glass
- Polyvinylchloride
- Plexiglas®
- Vinyl

AREAS OF USE:

Counters/Countertops, Operating Room Lights/Tables/Surfaces, **Dental** Instruments/Implements/Surfaces, **Medical** Instruments/Implements/Surfaces, **MadaJets, Ambulance** Equipment/Surfaces, **Laboratory** Equipment/Surfaces, **Respiratory Therapy** Equipment, **Respirators**/Face Pieces, **CPR Training** Mannequins, **Nail Care** Implements, **Exercise** Equipment/Benches/Mats, **Tanning Beds**, Hardhats, **Headphones**,

ORGANISMS	KILL TIME
Herpes Simplex II, HIV-I (Aids Virus), Influenza Strain A2/HK	30 sec
Rhinovirus, Trichophyton mentagrophytes	2 min
Methicillin Resistant Staphylococcus Aureus (MRSA)	3 min
Hepatitis B & C Virus (HBV, HBC), Pseudomonas aeruginosa, Salmonella enterica, Staphylococcus aureus	3 min
Aspergillus niger, Candida albicans, Vancomycin-resistant Enterococcus faecalis (VRE)	5 min
Mycobacterium bovis BCG (Tuberculosis)	6 min

ULTRASONIC CLEANING:

Thoroughly pre-rinse instruments under running water to remove visible gross debris. Using 1 ounce of this product per liter of water in ultrasonic unit, immerse instruments into solution and activate ultrasonic unit for 5 minutes (or longer if necessary). Remove and rinse instruments. Follow with appropriate cleaning and disinfection process. Change solution as needed when it becomes diluted or visibly soiled.