

Product Name: **Cultubes™**

Catalogue No.: **T406-2A**

This document replaces any previous version

1. Product Description:

- **14 ml Culture Tube with Cap:** Disposable, sterile , 17 x 95 mm, culture tube with graduations and white marking area

2. Packaging:

- **Case:** Packaged in 20 trays of 25 / 500 units per case

3. Product Specifications:

- **Material:**
 - Tube: Polypropylene
 - Cap: Polyethylene
- **Dual-position snap cap:**
 - Sealed position for anaerobic cultures
 - Loose position for aerobic work
- **Certified Pyrogen Free**
- **Gamma radiation sterilized** at a SAL of 10^{-3} ; specified dose between 6.5 kGy and 13.5 kGy
- **Temperature range:** -90°C to +100°C
- **Not autoclavable**
- **Centrifuge resistant** at up to 3 000 g

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resins conform to FDA 21 CFR 177.1520
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals
- **REACH (SVHC):** Plastics are in conformity to REACH standards
- **LATEX:** Materials are Latex Free
- **BSE / TSE:** Materials are BSE / TSE Free

5. Quality Assurance:

- **Clear,** no presence of contamination in plastic
- **Visual attributes**
- **Dimension and volume measurements**
- **Closure verification**

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination
- Protect from any damage to the packaging which could compromise the product sterility

8. Recommended Use:

- Verify proper cap closure when using biohazard material and / or chemical reagents
- Follow chemical resistance chart recommendations
- For use in automated equipment, follow the equipment manufacturer's instructions

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*

Product Name: **Cultubes™**

Catalogue No.: **T406-2**

This document replaces any previous version

1. Product Description:

- **14 ml Culture Tube with Cap:** Disposable, sterile 17 x 95 mm culture tube with graduations and white marking area

2. Packaging:

- **Case:** 20 trays of 25 units / 500 units per case

3. Product Specifications:

- **Material:**
 - Tube: Polystyrene
 - Cap: Polyethylene
- **Dual-position snap cap:**
 - Sealed position for anaerobic cultures
 - Loose position for aerobic work
- **Certified Pyrogen Free**
- **Gamma radiation sterilized** at a SAL of 10^{-3} ; specified dose between 6.5 kGy and 13.5 kGy
- **Temperature range for tube only:** -80°C to +80°C
- **Not autoclavable**
- **Centrifuge resistant** up to 1 400 g

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resins conform to FDA 21 CFR 177.1640 & CFR 177.1520
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals
- **REACH (SVHC):** Plastics are in conformity to REACH standards
- **LATEX:** Materials are Latex Free

5. Quality Assurance:

- **Clear, no presence of contamination** in plastic
- **Visual attributes**
- **Dimension and volume measurements**
- **Closure verification**

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any damage to the packaging which could compromise the product sterility

8. Recommended Use:

- Verify proper cap closure when using biohazard material and / or chemical reagents
- Follow chemical resistance chart recommendations
- For use in automated equipment, follow the equipment manufacturer's instructions

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*