



BINAXNOWTM STREP A

When the result is critical, so is the test.
Rapid in-office detection of Strep A
pharyngitis facilitates prompt and
targeted antibiotic treatment.




BinaxNOW™ Strep A

ABOUT

Detects *Streptococcus pyogenes* Group A antigen from throat swab specimens

- Simple, rapid test
- Highly accurate, sensitive and specific
- Easy to interpret
- Result in 5 minutes

DETECTION OF GROUP A STREP – BOTH SIMPLE AND ACCURATE



See the Product Insert for the full test procedure.

PERFORMANCE

BINAXNOW™ STREP A TEST VERSUS CULTURE

Sensitivity	92%
Specificity	100%
Overall Agreement	98%

ORDERING INFORMATION

PRODUCT	ORDER NO.
BinaxNow™ Strep A Test, 25 Test Kit	735-025
Control Swabs	730-010

OTHER BINAXNOW™ TESTS

BINAXNOW™ INFLUENZA A & B (not available in the United States)

BINAXNOW™ LEGIONELLA

BINAXNOW™ RSV

BINAXNOW™ S. PNEUMONIAE

BINAXNOW™ MALARIA

BINAXNOW™ G6PD

ABBOTT.COM/POCT

Abbott Rapid Dx International Limited
Parkmore East Business Park
Ballybrit, Galway
Ireland
Tel: +353 (0) 91 429 900

© 2020 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. Any photos displayed are for illustrative purposes only. Any person depicted in such photos is a model. 12000574-03 10/20

