

Introducing

Inject-Safe Barrier Bandage

The *first and only* bandage designed specifically for injections

- Improved Safety
- Lower Cost
- More Convenient
- Less Pain



Simple Procedure

Step 1: Prepare injection site with alcohol swab

Step 2: Remove Bandage from package and apply to skin

Step 3: Make injection through center of bandage remove and safely dispose of syringe



Improved Safety

- Barrier Bandage creates an antiseptic puncture site supporting 'no touch' technique
- Eliminates need to return to puncture site post injection
- Self sealing membrane contains bleeding and splattering during an injection
- Protects user from coming into contact with potential blood borne pathogens
- Simple procedure reduces potential for accidental finger stick injuries
- 98% of users stated that the Barrier Bandage made them feel safer

Lower Cost

- Combines the functionality of a band aid, cotton ball and rubber gloves
- More efficient procedure lowers overall cost of injection by as much as 50%

More Convenient

- Simple procedure saves time and integrates easily into busy office, clinic or pharmacy

Less Pain

- Unique concentric adhesion reduces injection pain for patients
- 70% of patients felt the Barrier Bandage made the injection less painful



Injection Safety Industries -1945 Scottsville Road B-2 #140 -Bowling Green, KY 42104

www.injectsafebandages.com

info@injectsafebandages.com