

Pro-Surface® Disinfectant Spray



TECHNICAL DATA SHEET

FEATURES & BENEFITS

- 1-minute surface contact time accelerates workflow and reduces patient transition time.
- Broad spectrum bactericidal activity eliminates 48 pathogens, including bacteria, viruses and fungi.
- Low toxicity profile will not damage non-porous instruments or surfaces; no special handling/protective eyewear required.
- This fragrance-free formulation increases compliance without disrupting office operations.
- Low quaternary ammonium content does not leave residue on surfaces.

RECOMMENDED USE

- Used to clean or disinfect hard, non-porous surfaces in dental offices, hospitals, and other healthcare settings.

PRODUCT DETAILS

Product Name:	Medicom Pro-Surface® Disinfectant Spray with TotalClean™ Technology
Product Type:	Disinfectant Spray
Active Ingredient:	Ethyl Alcohol 72.50 % w/w, Didecyl Dimethyl Ammonium Chloride 0.33% w/w, Other Ingredients 27.17 % w/w
Mililitres per Bottle:	946 ml
Country of Origin:	USA



CHARACTERISTIC	SPECIFICATION	
Material Composition	Bottle	High-density polyethylene (HDPE)
	Child Proof Cap	High-density polyethylene (HDPE)
	Seal	High-density polyethylene (HDPE)
	Label	Film
Description	Spray bottle containing disinfectant	
Color	Clear liquid	

CERTIFICATION & STANDARDS

ISO 13485 certified manufacturer.
Health Canada registered medical device manufacturer.
FDA Class EPA, EPA 88494-3-93710.

STORAGE CONDITIONS

Protect from extreme cold or heat.
Shelf life of 3 years.



Appendix 1 (List of items)

Item Code	Description	Packaging Information	Size	Color	Treatment	Marketed in USA as	Marketed in Canada as
40066	Medicom Pro-Surface Disinfectant Spray	12 Bottles/Case	32 fl. oz. (946 ml)	Clear Liquid	Non-Sterile	EPA 88494-3-93710	Not Marketed In Canada

AMD Medicom Inc. | 2555 Chemin de l'Aviation, Montreal, Quebec, Canada H9P 2Z2