


Characteristics

For maximum protection in high risk environments. Chlorinated and coated for rapid donning in critical situations. Textured finish with a medium surface tack for positive grip. High tensile strength for added protection. Tested for use with chemotherapy drugs. Blue color.



**Exam Glove
Non-Sterile**



**DermAssist®
EP BLUE®**

**Latex
Extra Protection
Series 181**

Extended Cuff

Features:

- Thicker and Longer than Regular Exam Gloves for Extra Protection in High Risk Situations
- Textured Finish for an Improved Wet/Dry Grip
- Low Protein for Reduced Chances of Latex Allergic Reactions
- Low Modulus for a Softer, More Comfortable Fit



PRODUCT DETAILS

SIZE	ITEM NO.	PACKAGING	DESCRIPTION
S	181100	50 Gloves/box, 10 boxes/case	Gloves, Exam, Latex, Chemo, Non Sterile, Powder-Free, Textured, Extended Cuff, Blue Color, 15 mil Finger Thickness
M	181200	50 Gloves/box, 10 boxes/case	
L	181300	50 Gloves/box, 10 boxes/case	
XL	181350	50 Gloves/box, 10 boxes/case	
XXL	181400	50 Gloves/box, 10 boxes/case	

View this product on
our website:



Product Solutions You Trust



INNOVATIVE
HEALTHCARE
CORPORATION

www.ihcsolutions.com

Specification Sheet

DERMASSIST[®] EP BLUE[®]

Powder-Free Latex Exam Gloves



- 50 Qty (By Weight)
- Single Use
- Ambidextrous
- Non-Sterile

DermAssist[®] EP Blue[®] is manufactured in compliance with multiple international standards, including the following:

Designation	Standard
ASTM D3578	Standard Specification for Rubber Examination Gloves
ASTM D5151	Standard Test Method for Detection of Holes in Medical Gloves
ASTM F1671	Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens
ASTM D5712	Standard Test Method for Analysis of Aqueous Extractable Protein in Natural Rubber

Average Length	Average Palm Thickness	Average Finger Thickness
11.5 in ± 290 mm	13.0 mil ± 0.33 mm	15.0 mil ± 0.38 mm

Tensile Strength & Elongation	Before Aging	After Accelerated Aging
Tensile Strength (Mpa)	28	24
ASTM Requirement Min. (Mpa)	18	14
Elongation (%)	850	820
ASTM Requirement Min. (%)	650	500

INNOVATIVE
HEALTHCARE
CORPORATION
Product Solutions You Trust
1-800-272-1533



Intertek

Innovative Healthcare Corporation
is certified to both ISO 13485:2003
QMS for medical devices.

Chemotherapy Drug Permeation (Breakthrough detection time in minutes, 0.01µg/cm ² /min.)	(ASTM D6978) Breakthrough Detection Time
Cisplatin (1.0 mg/mL)	>240
Cyclophosphamide (Cytoxan) (20.0 mg/mL)	>240
Dacarbazine (DTIC) (10.0 mg/mL)	>240
Doxorubicin Hydrochloride (2.0 mg/mL)	>240
Etoposide (20.0 mg/mL)	>240
Ifosfamide (50.0 mg/mL)	>240
5-Fluorouracil (50.0 mg/mL)	>240
Methotrexate (25.0 mg/mL)	>240
Mitomycin C (0.5 mg/mL)	>240
Mitoxantrone (2 mg/mL)	>240
Paclitaxel (Taxol) (6.0 mg/mL)	>240
Vincristine Sulfate (1.0 mg/mL)	>240
*Carmustine (BiCNU) (3.3 mg/mL)	DO NOT USE
*Thio-Tepa (10.0 mg/mL)	DO NOT USE

*Please call Customer Service at 1 800 272 1533 for a list of products appropriate for use with this drug.